

AmpliconBiSeq: User Guide

Altuna Akalin
altuna.akalin@fmi.ch

September 2, 2013

Contents

1	Introduction	1
2	Basics	1
2.1	Aligning reads	1
2.2	Conversion quality check by spike-ins	2
2.3	Extracting amplicon based methylation from alignments	3
2.4	Visualizing Amplicons	3
3	Calculation of meta-methylation profiles	4
4	Calculation of similarity between CpG pairs from the same amplicon	7
5	Convenience functions	7
5.1	Getting methylation ratio for all bases covered in experiments	7
5.2	Subsetting AmpliconViews	9
5.3	Getting amplicon and sample information	9
5.4	Other functions	9
6	Future Work	9

1 Introduction

AmpliconBiSeq is an R package for visualization and analysis of amplicon bisulfite sequencing experiments using high-throughput sequencing techniques. This is a targeted version of HT bisulfite sequencing. Amplicon bisulfite sequencing provides very deep coverage for short regions of the genome.

2 Basics

The entry point to an AmpliconBiSeq analysis is typically a set of aligned reads. These aligned reads will be processed by the package and an object that contains a summary for each amplicon of each sample will be produced. This basic object can be visualized with functions provided in the package and parts of it can be extracted for further analysis.

2.1 Aligning reads

Alignment of reads should be done by `QuasR` package. The entry point to AmpliconBiSeq analysis will be a `qProject` object obtained from `QuasR`. Below is an example code chunk that will achieve such alignment using `QuasR` package.

```

library(QuasR)
# make clusters, needed to run in parallel
cluObj=makeCluster(8)
#change this to desired location of your alignments
my.alignmentsDir="/work2/gschub/altuna/projects/AmpliconTimeCourse_cre_Feldmann/aln"
#assuming text files are in your working directory
proj=
  qAlign(
    sampleFile="~/w2/projects/AmpliconTimeCourse_cre_Feldmann/readFiles.v0.txt",
    genome="BSgenome.Mmusculus.UCSC.mm9",
    auxiliaryFile="/work2/gschub/Juliane/scripts/Altuna/auxFile.forQuasR.txt",
    aligner="Rbowtie",
    paired="fr",
    bisulfite="undir",
    projectName="AmpliconBiSeq",
    alignmentsDir=my.alignmentsDir,
    clObj=cluObj
  )

```

sampleFile and auxiliaryFile are tab separated text files that contain the sample names and read locations. For the auxiliaryFile, it contains locations of spike-in sequences that will also be aligned against. See ?qAlign for more information on sampleFile and auxiliaryFile requirements.

2.2 Conversion quality check by spike-ins

The qProject object that contains the aligned reads can also contain alignments for spike in sequences. These sequences will have predefined methylation properties and bisulfite conversion efficiency can be interrogated using those. The spikeCheckfunction provides methylation statistics for spike-ins. It plots a histogram or set of histograms from spike-in experiments which are helpful to deduce conversion efficiency of the experiment.

```

library(AmpliconBiSeq)

## Loading required package: grid
## No methods found in "Rsamtools" for requests: readBamGappedAlignments


spikeCheck(proj, auxName = "T7", sampleName = "mock4")

## Ratio of methylated Cs to Total number of Cs
##
## T7 : 0.9822
## $T7
## [1] 0.9822

```

T7 : mock4 Methylation Distribution

mean: 0.98 median: 1



2.3 Extracting amplicon based methylation from alignments


The `AmpliconViews` function extracts the amplicon related information from the alignment and calculates meta-methylation profiles and similarity between CpGs, as well as average methylation values and coverage per base. The function returns an `AmpliconViews` object. The object holds necessary information for further analysis and visualization.

```
ampliconLoc <- GRanges(seqnames = c("chr18", "chr18"),
  ranges = IRanges(c(69674375, 69674975), c(69674775,
    69675375)))
x <- AmpliconViews(proj, range = ampliconLoc, tag = c("amp1",
  "amp2"), sampleNames = c("mock4"), call.matrix = FALSE,
  verbose = TRUE)
x
```

2.4 Visualizing Amplicons

Amplicons can be visualized with `plotAmpliconView` function. The function can visualize one amplicon at a time, so it is important that users select which amplicon they want to visualize before hand. This can be done with `getAmplicon` function. An example of this is shown below.

```
data(ampViewEx) # load example data
myAmp <- getAmplicon(ampViewEx, "mock4", "chr18_69674375_69674775")
plotAmpliconView(myAmp)
```




3 Calculation of meta-methylation profiles

An amplicon can be represented as a binary methylation call matrix where each row is a read fragment and the columns are CpGs. Meta-methylation profiles can be thought of as patterns that explain a sub-section of methylation call matrix from an amplicon. This matrix can be analyzed to get methylation profiles that explains most of the individual profiles. This way, huge matrices can be efficiently visualized and summarized, and this could also be used as a metric for sample heterogeneity. Meta-methylation profiles are calculated using two rounds of singular-value decomposition (different way of doing PCA) and clustering. First round

of SVD removes the noise using a user defined cutoff or an 'automatic' cutoff. The cutoff designates what percentage of variation is considered as noise. The first round will basically remove principal components that does not explain much of the variation. The second round of SVD will be run on noise filtered matrix of methylation calls, and rows of the matrix will be clustered based on the top contributing components (the components that explain most variation), then for each cluster a meta-methylation profile is calculated by taking the average of methylation scores and binarizing them. The size of the cluster can be used to calculate what percentage of the data has driven the given meta-methylation profile.

In the example below, we will simulate a methylation call matrix which has three methylation profiles and some noise. Here are the methylation profiles:



Now, we will construct a simulated methylation call matrix and add some noise. "Profile1" will be replicated 500 times, "Profile2" will be replicated 300 times and "Profile3" will be replicated 200 times. Then we will add 5 % noise which will represent methylation call error or other noise.


**methylation call matrix
before noise addition**



**methylation call matrix
after noise addition**



Here are the returned meta-methylation profiles after analyzing the simulated call matrix with the method explained above. The method also orders the meta-methylation profiles based on how much of the data set they explain. In this example, "Meta1" returned as most important since it, explains 44% of the dataset. The method returned all the implanted profiles correctly and ordered them correctly based on how much of the data they can explain.



4 Calculation of similarity between CpG pairs from the same amplicon

For a given amplicon, a similarity between a pair of CpGs can be calculated as simple Jaccard similarity, which is number of common methylation calls over all methylation calls for a pair of CpGs.

$$JaccardSimilarity = \frac{M_1 \cap M_2}{M_1 \cup M_2} \quad (1)$$

In equation (1), $M_1 \cap M_2$ is the number of identical methylation calls for CpG_1 and CpG_2 , and $M_1 \cup M_2$ is the number of methylation calls that are present both CpGs. For example, the Jaccard similarity between these two binary vectors, $A=[1,1,1,0,0]$ and $B=[1,1,0,0,0]$ is 0.8 since 4 out of 5 elements are identical.

There are two other experimental similarity measures calculated. One is a method suggested by Landan et. al (doi:10.1038/ng.2442, Nature Genetics). The other one is still the Jaccard similarity but calculated after removing data points where both CpGs are unmethylated.

5 Convenience functions

These functions help query and extract parts of AmpliconViews object, so those parts can be utilized in further analysis.

5.1 Getting methylation ratio for all bases covered in experiments

`methRatio` function can be used to get a table of methylation ratios and coverage per base or per amplicon.

```
data(ampViewEx)
methRatio(ampViewEx)

## GRanges with 546 ranges and 5 metadata columns:
##      seqnames      ranges strand
```

```

##          <Rle>                <IRanges> <Rle>
##      [1] chr18 [25826977, 25826977]      *
##      [2] chr18 [25827148, 25827148]      *
##      [3] chr18 [69674494, 69674494]      *
##      [4] chr18 [69674628, 69674628]      *
##      [5] chr18 [69674644, 69674644]      *
##      ...
## [542] chr8 [114652776, 114652776]      *
## [543] chr8 [114652802, 114652802]      *
## [544] chr8 [114652823, 114652823]      *
## [545] chr8 [114652854, 114652854]      *
## [546] chr8 [114652905, 114652905]      *
##      | mock4.coverage mock4.methRatio
##      |      <numeric>      <numeric>
##      [1] |      3683      0.02851
##      [2] |      3694      0.03790
##      [3] |      9669      0.79150
##      [4] |      9772      0.95405
##      [5] |      9765      0.91091
##      ...
## [542] |      3165      0.010427
## [543] |     124834      0.095847
## [544] |     124744      0.133153
## [545] |     124405      0.078759
## [546] |     123559      0.001554
##      cre10.coverage cre10.methRatio
##      <numeric>      <numeric>
##      [1]      8789      0.01422
##      [2]      8799      0.02193
##      [3]      1776      0.64640
##      [4]      1798      0.60623
##      [5]      1786      0.50280
##      ...
## [542]      1541      0.041531
## [543]     107710      0.050515
## [544]     107666      0.047183
## [545]     107508      0.019226
## [546]     106984      0.004225
##      tag
##      <factor>
##      [1] Constitutive_LMR-REST
##      [2] Constitutive_LMR-REST
##      [3] FMR
##      [4] FMR
##      [5] FMR
##      ...
## [542] Constitutive_LMR-REST
## [543] Constitutive_LMR-REST
## [544] Constitutive_LMR-REST
## [545] Constitutive_LMR-REST
## [546] Constitutive_LMR-REST
## ---
## seqlengths:
##      chr1 chr10 chr11 chr12 ... chr8 chr9 chrX

```



```
##      NA      NA      NA      NA ...      NA      NA      NA
```

5.2 Subsetting AmpliconViews

`getAmplicon` is the function for subsetting an `AmpliconViews` object.

```
data(ampViewEx)
myAmp <- getAmplicon(ampViewEx, sampleNames = "mock4",
  ampliconNames = "chr18_69674375_69674775")
```

5.3 Getting amplicon and sample information

`getAmpliconNames`, `getSampleNames` are the functions that retrieve amplicon names and sample names. `getAmpliconRanges` returns `GRanges` object that contains the locations of the amplicons.

```
data(ampViewEx)
x <- getAmpliconNames(ampViewEx)
head(x)

## [1] "chr18_69674375_69674775"
## [2] "chr18_69674975_69675375"
## [3] "chr18_69675575_69675975"
## [4] "chr6_113705515_113705915"
## [5] "chr11_68768916_68769316"
## [6] "chr12_112491293_112491693"

getSampleNames(ampViewEx)

## [1] "mock4" "cre10"
```

5.4 Other functions

There are a number of other convenience functions to access particular parts of the `AmpliconViews` object. `getAvMeth` and `getCoverage` gets methylation per base and coverage per base for a specific amplicon. `getExampleMethMat` returns a sample of methylation call matrix (if there was enough coverage) on the amplicon. `getMethMat` returns full methylation call matrix for a given amplicon, if it was extracted during object creation.

6 Future Work

The package will support visualization of extra tracks `plotAmpliconView` function.